

<http://www.hwm.kr>



CPM (MC-2000)



Features

- ⇒ Special functions
 - Vibration : Provides muscle stimulation at end of ROM.
 - Progress : Increases ROM gradually until desired level achieved.
 - Limited Pause : Pauses up to 5 seconds at end of ROM.
- ⇒ Seven speeds are available for precise therapy options.
- ⇒ Hand and main control units are easy to use for patient and caregiver.
- ⇒ Timer can be set for up to 60 minutes.
- ⇒ Beeper switch notifies patient whether the button is clicking clearly.
- ⇒ Structure is convenient to carry, store and maintain.
- ⇒ Upper frame accommodates a wide range of patient sizes from children to large adults.
[Reference, 'Size Korea, Human, Engineering Data']
- ⇒ Handy control memory

Indications

- ⇒ Sports injuries
- ⇒ Total knee replacements
- ⇒ Open meniscectomies
- ⇒ ACL reconstructions
- ⇒ Joint debridement
- ⇒ Total hip replacements
- ⇒ Open reduction and internal fixation of intra-articular fractures.

Specifications

- ⇒ Range of motion $-5^{\circ} \sim 140^{\circ}$ deg
- ⇒ Adjustability
 - Femur size : 25 ~ 45 cm (9.8~17.7 inch)
 - Tibia size : 38 ~ 56 cm (14.9~22 inch)
- ⇒ Weight : 14 kgf (30.8lbf)
- ⇒ Length : 93 cm (36.6 inch)
- ⇒ width : 35 cm (13.8 inch)
- ⇒ Speed : Min (137° /min.) ~ Max (191° /min.) 7 STEP