

## METAMAX – The Mobile CPX System

Page: 2/3



## METASOFT – the metabolic Software



*click to enlarge*

### STRAIGHT-FORWARD OPERATION

- Customizable online instruction
- Configurable personal user profiles
- Data export in multiple formats (e.g. Excel, ASCII)
- AUTOMATED Tests
- Maximum oxygen uptake – VO<sub>2</sub> max

### Oxygen Kinetics

- Determination of the anaerobic threshold (AT) based on the V-Slope or Ventilatory Equivalent method

### CUSTOMIZED SCREENS AND REPORTS

- Variable adaptation of the parameters, test and analysis windows
- Print reports of your choice: summary report, UCLA (Wasserman) report, full disclosure formats - or individual the formats.

### UNPARALLELED Software UPGRADES

- Spirometry incl. flow volume loops during exercise
- Comprehensive Nutritional Assessment
- Cardiac Output
- CORTEX STATERA® - the ultimate tool for Patient Risk Stratification