


**Microwave traction bed:** Model Number: YHZ-2(microwave new type)



See larger image: Microwave traction bedcolour: green

Brand Name: xiangyu

Place of Origin: Henan China (Mainland)

Model Number: YHZ-2(microwave new type)

Port: Qingdao

Payment Terms: L/C, T/T

1. Microcomputer controlled lumbar traction + microwave physiotherapy combined treatment table, dual purpose for one machine.

1. Microcomputer controlled lumbar traction + microwave physiotherapy combined treatment table, dual purpose for one machine.

2. Traction therapy combines the microwave diathermy therapy, increase the treatment effect.

3. The lumbar diathermy can slide with the bed board during the traction. Increase the treatment effect greatly.

4. Intelligent digits display the time, traction force and the microwave power.

5. with the function of the traction force compensation.

6. There are continuous traction, intermission traction, repeat traction and other traction modes.

7. The lumbar microwave diathermy therapy (output power can be controlled by computer program, automatic adjustable microwave time.