

**NEW
MODEL**

BASISPLUS

40-MINUTES SESSION

**...FOR PREVENTION, ANTI-AGING, WELL-BEING,
SLEEP QUALITY AND RELAXATION.**


The new AIRNERGY Basis Plus is equipped with 2 AIRNERGY activation units and colour light therapy (blue). The recommended application time is 40 minutes.

EQUIPMENT	
2 AIRNERGY catalysts (activation units)	
Integrated colour light therapy (blue)	
3 intensity levels (50%, 75%, 100%)	
LCD display	
Five-language menu guidance (DE, EN, FR, IT, ES)	
Microcontroller-controlled electronics	
Individual programme settings	
Automatic self-test	
Switchable acoustic signal	
Optional: aromatherapy	

TECHNICAL SPECIFICATIONS	
Power supply	12 V DC
Power supply unit	100-240V AC/50-60Hz out: 1.25 A
Power consumption	max. 1,100 mA
Air supply	approx. 4 l per min.
Dimensions (mm) incl. glass bottle	350 x 310 x 240 (L x W x D)
Weight	3.1 kg

INFORMATION	
SERVICE	Every 2,000 hours of operation (approx. 3,000 sessions)
WARRANTY	2 years

SHIPMENT INCLUDES	1 power supply unit, 10 nasal cannulas, 1 glass bottle incl. sparkling element, 2 air filters, 1 cannula bag, 1 user manual
ACCESSORIES (OPTIONAL)	Aroma bottle, stainless steel support for aroma bottle
ARTICLE-NO.	2012055



2011_02_AIR-DB-PEN_AIR01

