

# Dental Chair SK 652 L



## Technical parameters:-

Power supply	230 V $\pm$ 10%
Frequency	50 Hz $\pm$ 5%
Total input	150 VA $\pm$ 10%
Lowest chair seat position	420 $\pm$ 20 mm
Highest chair seat position	850 $\pm$ 20 mm
Weight	95 KG $\pm$ 5 KG
Type of electrical protection	Class I.
Level of electrical protection	device type B
Environment:	as per the EN 60 601-1
Temperature range	
relative humidity range	
atmospheric pressure range	from +10 °C up to +40 °C
from 30 % up to 75 %	
from 700 hPa up to 1 060 hPa	