

PERCEPTION II



State-of-the-art technology at an affordable price.

Use the Perception II to monitor your lab area, exam rooms, or any climate-dependant area. Temperature, barometric pressure, and humidity are shown simultaneously on one easy-to-read display. With just the touch of a key, you can change the temperature reading from Fahrenheit to Centigrade, or the barometric pressure reading from inches of mercury to millimeters or millibars. Its sleek design makes it ideal for use on your desk, in the lab, or in the exam room. The Perception II has a backlit display for easy viewing and requires no installation. Just plug it in, set the time, and adjust the barometric pressure to your altitude.

Perception II Functions

- Inside temperature from 32° to 140° F (0° to 60° C)
- Inside humidity
- Barometric pressure (with stored value) in inches, millimeters or millibars

Accuracy

- Temperature: $\pm 1^{\circ}\text{F}$, $\pm 0.5^{\circ}\text{C}$
- Barometric Pressure: $\pm 0.05''$; $\pm 1.3\text{mm}$; $\pm 1.7\text{mb}$
- Humidity: $\pm 5\%$

Perception II Accessories

- AC-power adapter
- Reversible mounting base for use on desk, shelf or wall
- Detailed instruction booklet
- 1 year warranty

MedGraphics Order Information

Part Number

- 536159-001
- 830600-001
(with interface cable)

Order Contact Information

- Fax: 651-379-8333
- Email: groupsaleservices@medgraphics.com

Call your MedGraphics supply representative today at 1-800-655-1133



All specifications subject to change without notice. Products may vary from that illustrated. Please contact your Medical Graphics representative for latest information, pricing and product availability.

©Copyright, Medical Graphics Corporation, 2006. MedGraphics is a registered trademark of Medical Graphics Corporation.

Part #060031-001 REV B