


Product Range of Rehabilitation



Heat Therapy

Electric Wax Therapy



Model: XYL-10

Shortwave Diathermy



Megapulse Senior 265

Continuous and Pulse

RESISTANCE EXERCISE BAND



Sanctband Low Powder Reduced Protein latex resistance band is developed according to the production method used in powder free latex gloves

Cold Pack MODEL:00-3102

Laser Therapy Unit



- Mobile units has large storage and chilling capacity.
- Requires no pumping
- plugs into any 110v outlet

EMG biofeedback-MAXTENS 2000



- Fully adjustable
- Convenience
- Sturdy construction with an attractive appearance

Electromyography (EMG)



(EP) systems are recognized worldwide for their superior amplifiers,stimulators, and outstanding signal quality.

Multifunctional Traction Bed



HLA-200 is an optimal therapy system for release pain by a suitable wavelength of diode laser, without any damaging on skin and, it can be selected for each application respectively ; 660nm, 830nm and 950nm.

VIBRATOR: 01M7 Standard



Isokinetic System



MPS21 is an automated electronic robotic dynamometer used in the computerized assessment of muscle evaluation, prevention function & training.

Abdominal Flexion & Back Extension



Model No. JYZ-III

-Gone are clumsy weight stacks & primitive hydraulic Cylinders

-Features our smooth frictionless Magnatech resistance system

-No unnecessary machine configurations or adjustments to

confuse users or slow training

Gaitkeeper Rehab Treadmill



Improve the ergonomics of the therapist facilitation of an involved patient walking on the treadmill by raising the treadmill walking surface.

The Petite Basic System Model:701



-Sequential compression with four overlapping compression cells to ensure distal to proximal flow and comfortable compression

SAM'S INTERNATIONAL

Head Office:

BA-5A, ASHOK VIHAR, PHASE-I,

ASHOK VIHAR, DELHI-110052, INDIA

TEL:27213471, FAX: 011-27130045, MOB: 0-9971358977, 0-9811318477

E-MAIL: samsinternational@yahoo.com, ssehgal09@yahoo.com.sg

WEBSITE: www.sams-international.com

