

# Sanctband



## Resistive Exercise Band

To address concerns of skin irritation and discomfort caused by powder and protein in natural latex, Sanctband has developed a very Low Powder and Reduced Protein exercise band. Thus, minimizing discomfort during handling and minimizing the possibility of allergic reactions. Sanctband Resistive Exercise Band is an effective exercise product in physiotherapy rehabilitative and strength training.

Available in:

- a** 3 in 1 combo dispenser.
- b** 30 pcs dispenser.
- c** 6 yards roll dispenser.
- d** 50 yards roll dispenser.



Sanctband Resistive Band	Pull Force Test (Kg)	
	100%	300%
PEACH - Extra Light	Min 1.3	Min 2.5
ORANGE - Light	Min 2.1	Min 3.7
LIME GREEN - Medium	Min 2.6	Min 4.5
BLUEBERRY - Heavy	Min 3.2	Min 5.5
PLUM - Extra Heavy	Min 4.2	Min 7.7
GRAY - Super Heavy	Min 5.2	Min 9.2

## Resistive Exercise Tubing

Sanctband Resistive Exercise Tubing offers an alternative to the Resistive Exercise Band. Sanctband Resistive Exercise Tubing allows more variation of exercises, and is widely used both in the rehabilitative and strength training in various sports.

Available in:

- a** 25 ft Dispenser.
- b** 100 ft Dispenser.



Accessories Available :



Sanctband Resistive Tube	Pull Force Test (Kg)	
	100%	300%
PEACH - Extra Light	Min 2.1	Min 4.0
ORANGE - Light	Min 2.3	Min 4.6
LIME GREEN - Medium	Min 2.9	Min 5.8
BLUEBERRY - Heavy	Min 3.7	Min 6.9
PLUM - Extra Heavy	Min 4.2	Min 7.5
GRAY - Super Heavy	Min 4.8	Min 8.3

These Products Contain Natural Rubber Latex Which May Cause Allergy.

# Sanctband



# Sanctband

Professional Resistance Training  
Better in every way



[www.sanctband.com](http://www.sanctband.com)